



D I A L Y S I S P R O D U C T S

NIPRO | TP-SURE

TRANSDUCER PROTECTOR



RE-ORDER# TP-SURE+01

- Superior sterile pressure monitoring device for prevention of cross-contamination of the fluid pathway
- Inner filter is 0.1 microns, effective as a bacterial and viral barrier
- Individually packaged / 100 per box
- Individually capped for aseptic technique

SAFEGUARDING YOU... PROTECTING YOUR PATIENT



NIPRO MEDICAL CORPORATION
3150 N.W. 107th Avenue
Miami, FL 33172
Tel: (305) 599-7174
Fax: (305) 592-4621

For additional information on this product and/or the NIPRO CONNECT™ PROGRAM, call us toll-free at 1-800-647-7698/ext. 260
www.nipro.com

© 2008 Nipro Medical Corp. Printed in U.S.A.
TP-SURE™ is a trademark of Nipro Corporation.
November 2008

