




# DISPOSABLE DOSE-SAVER SYRINGE

## Minimize the waste. Maximize the savings.

- When compared with a competitive low dead space syringe, the Nipro 1cc Dose-Saver was shown to have 9% less "dead space" and the 3cc syringe 10% less dead space than the competitive syringe.\*
- Exclusive flat bevel on needle reduces aspiration of air bubbles.
- For precise volume measurement, highly translucent barrel features ml/cc graduations in increments of .25 and .5ml/cc, and elastomer stopper assures smooth, accurate motion along entire length of the barrel.



**NOW AVAILABLE**  
**1cc and 3cc**

\*Data on file, Nipro Medical Corporation.  
**Less dead space than competitors**  
*(blue indicates amount wasted after use)*

### Ordering Information

Permanently attached needle · Packed 100/box, 10 boxes/case

Order Number	Description	Order Number	Description
JD+01D2238	1cc w/22G x 1 1/2"	JD+03D2238	3cc w/22G x 1 1/2"
JD+01D2516	1cc TB Syringe w/25G x 5/8"		

Please consult product labels and inserts for indications, contraindications, hazards, warnings, cautions and directions for use.



# 10cc SYRINGE

## Safeguarding you. Protecting your patient.

- The **NIPRO** 10cc syringe is capped for improved infection control, to protect clinician and patient from accidental contamination.
- For greater convenience and procedural efficiency, Nipro 10cc syringe is compatible with most infusion pumps.
- Syringe offers all recognized Nipro quality features, including a clear, easy-to-read barrel and a specially formulated elastomer stopper for smooth, accurate motion along entire length of the barrel.



**CAPPED FOR SAFETY**

### Ordering Information

Graduation: 0.2cc (12cc Total) · Packed 100/box, 9 boxes/case

Also available in Luer Slip Tip with needle.

Order Number	Description	Order Number	Description
JD+10L-WEI	Luer Lock Tip	JD+10S-WEI	Luer Slip Tip
10cc LUER LOCK TIP W/NEEDLE			
JD+10L2025-WEI	20G x 1"	JD+10L2138-WEI	21G x 1 1/2"
JD+10L2038-WEI	20G x 1 1/2"	JD+10L2225-WEI	22G x 1"
JD+10L2125-WEI	21G x 1"	JD+10L2238-WEI	22G x 1 1/2"

Please consult product labels and inserts for indications, contraindications, hazards, warnings, cautions and directions for use.

