



Precision-crafted for delicate anatomy.



NIPRO SAFELET™ CATH

INTRAVENOUS ETFE CATHETER



Built for efficiency. Designed for comfort.

- All Safelet IV catheters are designed specifically to perform in delicate vasculature, delivering maximum efficiency and patient comfort.
- Smooth catheter-to-needle transition results in atraumatic insertion for optimum patient comfort.
- ETFE catheter material is biocompatible, minimizing risk of irritation and infection while catheter is indwelling.



Visibility for confident placement.

- Catheter is radiopaque and features a flashback "window" for quick visualization of blood.
- Hub is color-coded for rapid identification of catheter size.
- Available in a wide selection of sizes to match virtually any clinical need.



NIPRO SAFELET™ CATH

INTRAVENOUS ETFE CATHETER

User-friendly packaging with functional design.



- Hard plastic outer packaging protects the catheter before use.
- User-friendly construction of package prevents rolling on any work surface.

Order Information

Order#	Catheter	Hub Color	Order#	Catheter	Hub Color
CI+1450	14G x 2"	orange	CI+2032	20G x 1 1/4"	pink
CI+1632	16G x 1 1/4"	gray	* CI+2038	20G x 1 1/2"	pink
* CI+1650	16G x 2"	gray	* CI+2050	20G x 2"	pink
* CI+1663	16G x 2 1/2"	gray	CI+2225	22G x 1"	blue
CI+1832	18G x 1 1/4"	green	CI+2232	22G x 1 1/4"	blue
CI+1850	18G x 2"	green	CI+2420	24G x 3/4"	yellow
* CI+1863	18G x 2 1/2"	green	CI+2425	24G x 1"	yellow
CI+2025	20G x 1"	pink	CI+2450	24G x 2"	yellow

*Non-stock item.

Please consult product labels and inserts for indications, contraindications, hazards, warnings, cautions and directions for use.



NIPRO MEDICAL CORPORATION
 3150 N.W. 107 Avenue,
 Miami, FL 33172
 (305) 599-7174
 Fax: (304) 599-8454

© 2004, Nipro Medical Corp. Printed in U.S.A.
 SAFELET is a trademark of Nipro Medical Corp.
 0000 APA 0/00 OM