

NIPRO

Dialysis System

NIPRO SAFETOUCH FISTULA NEEDLE

Safety feature reduces risk of accidental needle puncture

Permanent safety lock

Pre-removal activation

Color coded for easy identification

Latex Free

Available in single and twin packs


Gauge Size

14G

15G

16G


17G



Engineered from the beginning to be a safety needle

NIPRO SAFETOUCH SAFETY FISTULA NEEDLE

FEATURES



Safetouch safety fistula needle

- Latex Free
- Safety feature reduces risk of accidental needle puncture
- Sterile
- Non-toxic, Non-pyrogenic
- Ultra thin needle walls allow maximum blood flow
- Ultra sharp for easier penetration
- Easy to use occlusion clamp
- Siliconized for smoother, more comfortable penetration
- Color-coded wings for easy identification
- Easy removal of Luer cap
- Luer tapered for easier connection with all blood lines
- Available with back eye for high flow rate, in fixed and rotating wing types
- Packaging: 50pcs/box, 10 boxes/case
- Available in twin pack (2 needles per package)

SPECIFICATIONS

REORDER CODE	DESCRIPTION	REORDER CODE	DESCRIPTION
FS+142530BC	14Gx 1"	FS+142530TP	14Gx 1"
FS+143230BC	14Gx 1 1/4"	FS+143230TP	14Gx 1 1/4"
FS+152530BC	15Gx 1"	FS+152530TP	15Gx 1"
FS+153230BC	15Gx 1 1/4"	FS+153230TP	15Gx 1 1/4"
FS+162530BC	16Gx 1"	FS+162530TP	16Gx 1"
FS+163230BC	16Gx 1 1/4"	FS+163230TP	16Gx 1 1/4"
FS+172530BC	17Gx 1"	FS+172530TP	17Gx 1"
FS+173230BC	17Gx 1 1/4"	FS+173230TP	17Gx 1 1/4"

The Above Reorder Codes Represents 12" Tubing, Fixed Wing and Backeye

All Sizes are Available With: Rotating Wing, 6 or 12 Inch Tubing, With or Without Backeye

Code #'s that End With "TP" Designates "Twin Pack" 2 Needles in Package

Distributed by:



NIPRO MEDICAL CORPORATION

3150 N.W. 107 Avenue
Miami, Florida 33172
Tel: (305) 599-7174
Fax: (305) 599-8454

Factories: Nipro Corporation, Ohdate Factory (Japan)
Nipro Medical Industries, Ltd. (Japan)
Nipro Corporation, Ltd. (Thailand)
Fuzhou Nipro Corporation, Ltd. (China)
Nipro (Shanghai) Co., Ltd. (China)
Nipro Medical LTDA (Brazil)

Overseas Office: N.V. Nipro Europe S.A. (Belgium, Holland, Milan, Madrid)
Nipro Medical Corporation (Miami, Argentina, Chile)
Nipro Medical Panama, S.A. (Panama)
Nipro Medical De Mexico, S.A. De C.V. (Mexico)
Nipro Medical Ltda. Sao Paulo Office (Brasil)
Nipro Corporation, Ltd. (Bangkok-Thailand)
Nipro Asia Pte. Ltd (Singapore, New Delhi, Dubai)
Nipro Corporation, Shanghai Liason Office (Shanghai-China)